

Eppendorf Certificate

Certificate of Purity - sterile

This package contains a high-quality consumable manufactured under the sterile Eppendorf Purity Standard.

The Eppendorf sterile consumables are produced in a controlled environment according to ISO class 8 of ISO 14644-1.

For this product Eppendorf certifies the following:

- > sterile
- > pyrogen-free



Quality control and subsequent certification is done by an independent laboratory accredited to ISO 17025. Eppendorf guarantees the conformity within the following limits:

Sterility	in accordance with USP, Ph. Eur. 2.6.12
Pyrogens (Endotoxins)	< 0.001 EU/mL (tested according to Ph. Eur. 2.6.14 (LAL test))

Lot-specific certificates are available on request or on the internet at www.eppendorf.com/lot-certificates.

Hamburg, April 2024

Joana Tziolis
Product Life Cycle Manager
Division Consumables

Thomas Hengstmann
Head of Global Quality Operations

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified

Eppendorf® and the Eppendorf Brand Design are registered trademarks of Eppendorf SE, Germany. All rights reserved, incl. graphics and images. Copyright ©2024 by Eppendorf SE.