



# Certificate of Containment Testing

**Containment Testing of Rotor FA-24x2 with  
spring-loaded lid construction (including a  
shortened seal)**

**in an Eppendorf 5425 Bench Top Centrifuge  
for use in Bench Top Centrifuges 5425 /R**

**Report No. 19/081**

**Report Prepared For:** Eppendorf AG, Hamburg, Germany

**Issue Date:** 15 April 2020

## Test Summary

Rotor FA-24x2 (5495HL303272) and with spring-loaded lid construction (including a shortened seal) was containment tested in an Eppendorf 5425/R bench top centrifuge, using Annex AA of IEC 61010-2-020:2016 (3<sup>rd</sup> Ed.). The sealed rotor was shown to contain a spill.

### Report Written By

**Name:** Ms Helen Hookway

**Title:** Biosafety Scientist

### Report Authorised By

**Name:** Mrs Sara Speight

**Title:** Senior Biosafety Scientist