


# Eppendorf Certificate

## Certificate of Purity - Biopur®

This package contains a high-quality consumable manufactured under the Biopur Eppendorf Purity Standard. The Eppendorf Biopur consumables are produced in an ISO class 8 clean room environment according to ISO 14644-1.

For this product Eppendorf certifies the following:

- > Sterile
- > Pyrogen/Endotoxin-free
- > Human DNA-free
- > Bacterial DNA-free
- > DNase-free
- > RNase-free
- > PCR inhibitor-free
- > ATP-free



These parameters are continuously monitored by an independent certified laboratory. Eppendorf guarantees the conformity within the following limits:

Sterility	in accordance with USP, Ph. Eur. 2.6.12
Pyrogens	< 0.001 EU/mL (tested according to Ph. Eur. 2.6.14 (LAL test))
Human DNA	< 2 pg
Bacterial DNA	< 50 fg
DNase	< 1.0 x 10 <sup>-6</sup> Kunitz units
RNase	< 1.0 x 10 <sup>-9</sup> Kunitz units
PCR inhibition	less than 10 targets amplifiable
ATP	< 1.0 x 10 <sup>-12</sup> mg

Lot-specific certificates can be downloaded from the internet at [www.eppendorf.com/certificates](http://www.eppendorf.com/certificates).

October, 2019

Joana Tziolis  
Expert Liquid Handling  
LH-Consumables

Dr. Birgit Schreiber  
Vice President Quality Management &  
Regulatory Affairs

**Your local distributor:** [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany  
E-mail: [eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com)

ISO 9001  
Certified

ISO 13485  
Certified

ISO 14001  
Certified