

Eppendorf Certificate

Certificate of Quality – ep Dualfilter T.I.P.S.®

ep Dualfilter T.I.P.S.® filter tips do neither differ in handling nor in specification for any of the properties listed below.

- > Orifice diameter, inner diameter, outer diameter
- > Length and geometry of the pipette tips
- > Weight / Attachment forces / Ejection forces
- > Compatibility with Eppendorf pipettes and non-Eppendorf pipettes
- > Quality of purity and sterility
- > Pipette tip material / Resistivity to chemicals
- > Position of the filter within the tip
- > Packing quantity, packing volume, packing weight

Volume	Order no.		
	old (international)	old (North America)	new
ep Dualfilter T.I.P.S.® Racks PCR clean and sterile			
0.1 - 10 µL S, 34 mm	0030 077.504	022491202	0030 078 500
0.1 - 10 µL M, 40 mm	0030 077.512	022491211	0030 078 519
0.5 - 20 µL L, 46 mm	0030 077.520	022491229	0030 078 527
2 - 20 µL, 53 mm	0030 077.539	022491270	0030 078 535
2 - 100 µL, 53 mm	0030 077.547	022491237	0030 078 543
2 - 200 µL, 55 mm	0030 077.555	022491296	0030 078 551
20 - 300 µL, 55 mm	0030 077.563	022491245	0030 078 560
50 - 1000 µL, 76 mm	0030 077.571	022491253	0030 078 578
50 - 1250 µL L, 103 mm	0030 077.750	022494002	0030 078 594
0.25 - 2500 µL, 115 mm	-	022491881	0030 078 586
0.2 - 5000 µL L, 175 mm	0030 077.725	0030077725	0030 078 624
100 - 5000 µL, 120 mm	0030 077.580	022491261	0030 078 616
0.5 - 10 mL L, 243 mm	0030 077.598	022491288	-
ep Dualfilter T.I.P.S.® Racks Forensic Grade			
0.1 - 10 µL M, 40 mm	0030 077.768	0030077768	0030 078 810
2 - 20 µL, 53 mm	0030 077.776	0030077776	0030 078 829
2 - 200 µL, 55 mm	0030 077.784	0030077784	0030 078 837
50 - 1000 µL, 76 mm	0030 077.792	0030077792	0030 078 845
ep Dualfilter T.I.P.S.® LoRetention® Racks PCR clean and sterile			
0.1 - 10 µL S, 34 mm	0030 077.610	022493000	0030 078 632
0.5 - 20 µL L, 46 mm	0030 077.628	022493002	0030 078 640
2 - 100 µL, 53 mm	0030 077.644	022493006	0030 078 659
2 - 200 µL, 55 mm	0030 077.741	0030077741	0030 078 667
20 - 300 µL, 55 mm	0030 077.636	022493004	0030 078 675
50 - 1000 µL, 76 mm	0030 077.652	022493008	0030 078 683
ep Dualfilter T.I.P.S.® 384 PCR clean and sterile			
0.1 - 20 µL, 42 mm	0030 078.403	0030078403	0030 078 853
0.5 - 100 µL, 52.8 mm	0030 078.411	0030078411	0030 078 861

Eppendorf Certificate

Volume	Order no.		
	old (international)	old (North America)	new
ep Dualfilter T.I.P.S.® SealMax Racks PCR clean and sterile			
0.1 - 10 µL S, 34 mm	0030 077.806	0030078691	0030 078 691
0.5 - 20 µL L, 46 mm	0030 077.814	0030078705	0030 078 705
2 - 100 µL, 53 mm	0030 077.822	0030078713	0030 078 713
2 - 200 µL, 55 mm	0030 077.830	0030078721	0030 078 721
20 - 300 µL, 55 mm	0030 077.849	0030078730	0030 078 730
50 - 1000 µL, 76 mm	0030 077.857	0030078748	0030 078 748

Hamburg, April 2022

Page 2 of 2



Joana Tziolis
Product Life Cycle Manager
Division Consumables



Monika Schneider
Vice President Global Quality Management &
Regulatory Affairs

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified

Eppendorf® and the Eppendorf Brand Design, ep Dualfilter T.I.P.S.® and LoRetention® are registered trademarks of Eppendorf SE, Germany. All rights reserved, incl. graphics and images. Copyright ©2022 by Eppendorf SE.