

Eppendorf Certificate

Certificate of Purity - sterile

This package contains a high-quality consumable manufactured under the sterile Eppendorf Purity Standard.

The Eppendorf sterile consumables are produced in an ISO class 8 clean room environment according to ISO 14644-1.

For this product Eppendorf certifies the following:

- > sterile
- > pyrogen-free



Quality control and subsequent certification is done by an independent laboratory accredited to ISO 17025. Eppendorf guarantees the conformity within the following limits:

Sterility	in accordance with USP, Ph. Eur. 2.6.12
Pyrogens (Endotoxins)	< 0.001 EU/mL (tested according to Ph. Eur. 2.6.14 (LAL test))

Lot-specific certificates are available on request or on the internet at www.eppendorf.com/certificates.

August, 2019

Joana Tziolis
Expert Liquid Handling
LH-Consumables

Dr. Birgit Schreiber
Vice President Quality Management &
Regulatory Affairs

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified