

## Press Release

# The New Eppendorf CryoCube® F740 Series

Hamburg, May 2017

Based on the reliable and durable Innova U725/ -G models, the new CryoCube F740 series by Eppendorf makes ULT sample storage future-proof.

An increased capacity of 740 L/ 576x boxes, improved energy efficiency without sacrificing temperature accuracy, and a considerable noise reduction create a safe storage environment. Options such as air-cooling or water-cooling, right-handed or left-handed doors, and three or five interior compartments set new benchmarks for user demands.

The CryoCube F740 is network-compatible and can be monitored through smart Eppendorf VisioNize® systems – directly at the device’s touchscreen or flexibly at a computer or tablet in the office.

### About Eppendorf AG:

Eppendorf is a leading life science company that develops and sells instruments, consumables, and services for liquid handling, sample handling, and cell handling in laboratories worldwide. Its product range includes pipettes and automated pipetting systems, dispensers, centrifuges, mixers, spectrometers, and DNA amplification equipment as well as ultra-low temperature freezers, fermentors, bioreactors, CO2 incubators, shakers, and cell manipulation systems. Consumables such as pipette tips, test tubes, microliter plates, and single-use bioreactor vessels complement the range of highest-quality premium products.

Eppendorf products are most broadly used in academic and commercial research laboratories, e.g., in companies from the pharmaceutical and biotechnological as well as the chemical and food industries. They are also aimed at clinical and environmental analysis laboratories, forensics, and at industrial laboratories performing process analysis, production, and quality assurance.

Eppendorf was founded in Hamburg, Germany in 1945 and has more than 3,000 employees worldwide. The company has subsidiaries in 26 countries and is represented in all other markets by distributors.

#### Press Contact:

Bettina Doerler  
Barkhausenweg 1  
22339 Hamburg

Tel: +49 40 538 01-640  
Fax: +49 40 538 01-662  
doerler.b@eppendorf.de

[www.eppendorf.com](http://www.eppendorf.com)